


BVI



Malosa® Single-Use Instruments

High Quality | Consistent | Efficient | Cost-Effective

bvimedical.com

Malosa[®] Single-Use Instruments

- **Eliminate Reprocessing and Sterilisation Costs**


No added burden on facilities to re-sterilise or repair damaged instruments, helping to increase operational efficiencies

- **Improve Patient Safety and Decrease Contamination**

By utilising a brand new instrument for each surgery, a patient's risk for infection or cross contamination is decreased

- **Consistency throughout all procedures**

Malosa[®] single-use instruments provide consistent and reliable performance for every case



C. Single-use instruments deserve consideration for eye surgery. Ocular Surgery News Europe Edition. December 2016

Southworth, P.M. Infections and exposures: reported incidents associated with unsuccessful decontamination of reusable surgical instruments. The Journal of Hospital Infection. September 17, 2014

Tarantola MD, R.; Shuler MD, PhD, M.; Kaushal, MD PhD, S.; Nagpal, MS, DO, FRCS(UK), M.; Disposable Instruments are Cost Effective and as Good as Reusable. Retina Today. Jan/Feb 2015

Table of Contents

Forceps	
Capsulorhexis Forceps	4
Fixation Forceps	5
Tying Forceps	6
IOL Forceps	6
General Forceps and Fixation Rings	7
Haemostatic Forceps	8
Dressing Forceps	8
Needle-Holders	
Corneal Needle-Holders	9
General Needle-Holders	9
Choppers	10
Double-Ended Choppers	11
Manipulators	12
Scissors	
Sprung Scissors	14
Dissection Scissors	15
Dressing Scissors	15
Speculums	
Adjustable Speculums	16
Sprung Wire Speculums	17
Paediatric and Neonatal Speculums	17
Hooks and Retractors	18
Depressors	18
Lacrimal	19
Markers and Gauges	20
Clamps	21
Curettes	21
Flap-Lifters and Corneal Markers	22
SMILE	23
Spatulas and Wells	23
Keratoplasty	24
Corneal Ring	24
Sclerectomy and Trabeculectomy	25
I-Ring	25
Miscellaneous	25



Note: Single-use instruments that are available with Contour Clips will be marked with a "CS" at the end of the product SKU.

Capsulorhexis Forceps

The Malosa range of single-use Capsulorhexis Forceps offers tip and blade designs for performing continuous curvelinear capsulotomy and supporting sub - 1.8 mm microincision (MICS) procedures.



Capsulorhexis Forceps	Description	SKU #	SKU # CLIPS
MICS Capsulorhexis - Curved/Sharp	0.4 mm x 0.3 mm Sharp tips, 12 mm Curved direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1567S	MMSU1567CS
MICS Capsulorhexis - Straight/Sharp	0.4 mm x 0.3 mm Sharp tips, 12 mm Straight direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1568S	MMSU1568CS
MICS Capsulorhexis - Straight/Grasping	0.4 mm x 0.3 mm Blunt tips, 12 mm Straight direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1571S	MMSU1571CS
MICS Capsulorhexis - Curved/Grasping	0.4 mm x 0.3 mm Blunt tips, 12 mm Curved direct action jaws, 1.2 mm at Base of shaft, 80 mm Flat handle	MMSU1573S	MMSU1573CS
Next Generation Capsulorhexis Forceps Sharp	Sharp curved tips, Round handle	581488	-
Next Generation Capsulorhexis Forceps Grasping	Grasping curved tips, Round handle	581486	-



MICS Capsulorhexis - Curved/Grasping



MICS Capsulorhexis - Straight/Sharp



Next Generation Capsulorhexis Forceps Grasping



Fixation Forceps



Fine Toothed Forceps Short (Bonn)




Notched Forceps (Hoskins)



Toothed Fixation Forceps	Description	SKU #	SKU # CLIPS
Fine Toothed Forceps Short (Bonn)	0.3 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1123S	MMSU1123CS
Toothed Forceps (Bonn)	0.45 mm 1 x 2 Teeth, 5 mm Tying platforms, 80 mm Handle	MMSU1124S	MMSU1124CS
Colibri Toothed Forceps Short	0.4 mm 1 x 2 Teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1129S	MMSU1129CS
Toothed Forceps 15 mm Shaft	0.4 mm 1 x 2 Teeth, 15 mm Shaft with 7.2 mm tying platforms, 80 mm Handle	MMSU1168S	MMSU1168CS
45° Toothed Forceps	0.3 mm 1 x 2 Teeth, 5 mm Tying platforms, 80 mm Handle	MMSU1139S	MMSU1139CS
45° Fine Toothed Forceps Short	0.27 mm 1 x 2 Teeth, 45 Degree teeth, 5 mm Tying platforms, 52 mm Handle	MMSU1145S	MMSU1145CS

Notched Fixation Forceps	Description	SKU #	SKU # CLIPS
Fine Notched Forceps Short (Hoskins)	0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle	MMSU1114S	MMSU1114CS
Notched Forceps (Hoskins)	0.3 mm Notched tips, 5 mm Tying platforms, 80 mm Handle	MMSU1115S	MMSU1115CS
Colibri Notched Forceps Short	0.2 mm Notched tips, 5 mm Tying platforms, 52 mm Handle	MMSU1209S	MMSU1209CS



Colibri Toothed Forceps Short

Tying Forceps



Tying Forceps	Description	SKU #	SKU # CLIPS
Straight Tying Forceps	0.35 mm Tips, 5 mm Tying platforms, 80 mm Handle	MMSU1125S	MMSU1125CS
Straight Tying Forceps Short	0.26 mm Tips, 5 mm Tying platforms, 52 mm Handle	MMSU1219S	MMSU1219CS
Curved Tying Forceps	0.38 mm Tips, 10 mm Curved tying platforms, 80 mm Handle	MMSU1210S	MMSU1210CS
Angled Tying Forceps (McPherson)	0.35 mm Tip, 5 mm Tying platforms, 8 mm Shaft angled 45 degree, 80 mm Handle	MMSU1127S	MMSU1127CS
Birks Tying Forceps	0.1 mm x 0.3 mm Tips, 5 mm Tying platforms, 80 mm Handle	-	MMSU1414CS
Waqar Suture Removal Forceps	0.3 mm x 0.65 mm Pick angled 90 degree, 5 mm Tying platforms, Curved 28 mm shaft, 80 mm Handle	MMSU1583S	-

IOL Forceps



IOL Forceps	Description	SKU #	SKU # CLIPS
IOL Implantation Forceps	7 mm Highly polished curved convex jaws, 80 mm Handle	MMSU1126S	MMSU1126CS
Lens Folding Forceps Rounded	Rounded inner jaws, 9.5 mm 45 Degree angled shaft, 80 mm Handle	MMSU1138S	-
Lens Folding Forceps Triangular	Triangular shaped inner jaws, 10.5 mm 45 Degree angled shaft, 80 mm Handle	MMSU1131S	-
Lens Loading Forceps	0.6 mm x 8 mm Tips, 45 Degree angled shaft, 80 mm Handle	MMSU1146S	MMSU1146CS
IOL Removal Forceps	0.1 mm x 0.1 mm Ridges along 4 mm of 15 mm shaft, Angled at 30 degree, 52 mm Handle	MMSU1460S	MMSU1460CS

General Forceps and Fixation Rings



General Forceps	Description	SKU #
Moorfields Forceps	0.6 mm tips, 11.5 mm Serrated jaws, 105 mm Long	MMSU0277S
Iris Pointed Forceps	0.5 mm × 0.5 mm Tips, 12 mm Serrated jaws, 105 mm Long	MMSU0247S
Iris Fine Toothed Forceps	1 mm × 1.15 mm 1 × 2 Teeth, 70 mm Handle	MMSU0278S
Fine Stem Toothed Forceps	1.65 mm 1 × 2 Teeth, 70 mm Handle	MMSU0080S
Max Fine Micro-Corneal Forceps	0.2 mm × 0.2 mm × 11 mm Curved tips, 70 mm Handle	MMSU1445S
Bishops Forceps	0.7 mm × 0.9 mm 1 × 2 Teeth, 70 mm Handle	MMSU0276S
St Martins Forceps	0.9 mm x 0.26 mm 1 × 2 Teeth, 52 mm Handle	MMSU0279S
Epilation Forceps	5 mm x 7 mm Flat plane tips, 70 mm Handle	MMSU0085S
Jewellers Forceps Fine Straight	7 mm Pointed tips, 80 mm Handle	MMSU0286S
Jewellers Forceps Fine Angled	7 mm, 45 Degree angled pointed tips, 80 mm Handle	MMSU0283S
Jewellers Forceps Angled/Smooth	4 mm 45 Degree angled pointed tips, 80 mm Handle	MMSU1472S
Jewellers Forceps Curved/Plain	Fine point "Swan Neck" curved shaft, 85 mm Handle	MMSU0284S
Scleral Plug Forceps	0.6 mm × 0.4 mm Tips, 9 mm Curved platforms, 70 mm Handle	MMSU1411S
Contact Lens Forceps	1.15 mm × 2.3 mm Tips with 0.5 mm notch, 6 mm 55 Degree angled shaft, 65 mm Handle	MMSU1284S
Adson Forceps, Toothed	1.35 mm × 1.8 mm 1 × 2 Teeth, 122 mm Long	MMSU0238S
Adson Forceps, Non-Toothed	0.85 mm × 1 mm Tips, 14.5 mm Serrated jaws, 122 mm Long	MMSU0241S

Fixation Rings

Fixation Ring	Description	SKU #
Fixation Ring	0.5 mm Thick 13 mm internal diameter ring, 1.5 mm × 0.6 mm Teeth at 3.5 mm intervals, 1.3 mm × 0.4 mm Teeth at 0.9 mm intervals, 75 mm Handle	MMSU1495S
Thornton Fixation Ring	0.5 mm Thick 6.5 mm radius ring, 2 mm × 0.9 mm Teeth at 5 mm intervals, 75 mm Handle	MMSU1137S



Fixation Ring

Haemostatic Forceps



Haemostatic Forceps	Description	SKU #
Halstead Mosquito Straight	2 mm Tip, 22 mm Serrated jaws without groove, 130 mm Long	MMSU0246S
Halstead Mosquito Fine Curved	1.1 mm x 0.85 mm Tips, 21 mm Serrated jaws, 130 mm Long	MMSU3007S
Halstead Mosquito Fine Straight	1.8 mm Tip, 22 mm Serrated jaws, 130 mm Long	MMSU3008S
Micro Cross-Action Bulldog Clamp	2 mm x 13.5 mm Grooved jaws, 42 mm Long	MMSU0199S

Dressing Forceps



Dressing Forceps	Description	SKU #
Straight Rampley	11 mm x 23 mm Serrated tips, 180 mm Long, Nylon plastic	MMSU0068S
Rampley Short	17 mm x 24 mm Serrated tips, 125 mm Long, Nylon plastic	MMSU1578S

Needle-Holders



Corneal Needle-Holder



Barraquer Needle-Holder



Troutman Needle-Holder



Corneal Needle-Holders

Corneal Needle-Holders	Description	SKU #	SKU # CLIPS
Corneal Needle-Holder	9 mm Fine curved jaws, 0.45 mm Tips, 106 mm Long	MMSU1110S	MMSU1110CS
Barraquer Needle-Holder	0.42 mm x 0.45 mm Tips, 4 mm Fine curved jaws with 5 mm shaft, 138 mm Long	MMSU3015S	-
Corneal Needle-Holder Long	0.45 mm Tips, 9 mm Curved jaws, Round handle, 138 mm Long	MMSU3017S	-

General Needle-Holders

General Needle-Holders	Description	SKU #
Troutman Needle-Holder	4 mm Fine curved jaws, 0.55 mm Tips, 76 mm Long handle	MMSU1243S
Castroviejo Needle-Holder Straight	1.5 mm Tips, 12.5 mm Straight jaws, Non-Locking, 142 mm Long	MMSU3009S
Castroviejo Needle-Holder Curved	1.5 mm Tips, 12.5 mm Curved jaws, Non-Locking, 142 mm Long	MMSU3011S
Castroviejo Needle-Holder Straight with Lock	1.5 mm Tips, 12.5 mm Straight jaws with lock, 142 mm Long	MMSU3012S
Castroviejo Needle-Holder Curved with Lock	1.5 mm Tips, 12.5 mm Curved jaws with lock, 142 mm Long	MMSU3013S



Troutman Needle-Holder

Choppers

The Malosa range of single-use choppers offers a range and combination of second instruments to suit all techniques of phacoemulsification.



Sadiq Chopper
MMSU1104S



Claw Chopper
MMSU1105S



Rosen Chopper
MMSU1106S



Phaco Chopper
MMSU1108S



Koch Chopper 1.9 mm
MMSU1118S




Koch Chopper 1.4 mm
MMSU1228S



**Koch Chopper 1.4 mm
Thin Shaft**
MMSU1440S




Nagahara Chopper
MMSU1441S



Denis Nucleus Chopper
MMSU1471S



**Greens Chopper
Right Handed**
MMSU1119S



**Greens Chopper
Left Handed**
MMSU1120S




Neuhann Nucleus Divider
MMSU1489S


Double-Ended Choppers




Phaco Chopper |
Mushroom Manipulator
MMSU1221S



Horizontal Chopper |
Kuglen Hook
MMSU1231S



Phaco Chopper |
Mushroom Manipulator Fine
MMSU1487S




Ursell Chopper | Manipulator
MMSU1238S




Diablo Chopper |
Flat Button Manipulator
MMSU1239S



Koch Chopper 1.4 mm |
Button Manipulator
MMSU1431S




Ayliffe Chopper | Chopper
MMSU1452S



Phaco Chopper |
Govan Manipulator
MMSU1496S

Manipulators

Malosa's portfolio of single-use manipulators offers a range and combination of second instruments to suit all surgical techniques.



Govan Manipulator
MMSU1102S



Drysdale Manipulator
MMSU1144S



Sinskey Hook 0.2 mm
MMSU1103S



Sinskey Hook 0.1 mm
MMSU1173S



Sinskey Hook Straight
MMSU1147S



Hagège Manipulator
MMSU1436S



Spatula
MMSU1291S



Button Manipulator
MMSU1142S



Mushroom Manipulator
MMSU1179S



Diablo Manipulator
MMSU1428S




Barraquer Iris Spatula
MMSU1457S



Butcher Rotator
MMSU1434S



Flat Round Tip Manipulator
MMSU1432S



Horizontal Y Manipulator
MMSU1204S



Vertical Y Manipulator
MMSU1109S



Lebuisson Manipulator
MMSU1458S



Lester Hook
MMSU1207S



Kuglen Hook Angled
MMSU1433S



**Kuglen Hook
Straight**
MMSU1225S



McPherson Manipulator
MMSU1143S



**Barraquer Iris
Spatula Round**
MMSU1592S



ICL Manipulator
MMSU1153S




0.1 mm Reverse Sinskey
MMSU1473S



Reverse Dialling Hook
MMSU1480S



Micro Manipulator
581441



Sinskey Hook
581442




Scissors

Sprung and direct-action ophthalmic scissors for all surgical specialties.



Westcott Scissors



Vannas Scissors



IOL Cutter



Castroviejo Scissors



Sprung Scissors

Sprung Scissors	Description	SKU #	SKU # CLIPS
Westcott Scissors Curved Pointed	0.25 mm x 0.3 mm Pointed tips, 12 mm Curved blades, 116 mm Long	MMSU1211S	-
Westcott Scissors Straight Rounded	0.5 mm x 0.6 mm Rounded tips, 17 mm Straight blades, 120 mm Long	MMSU1212S	MMSU1212CS
Westcott Scissors Curved Dolphin Nose	0.5 mm x 0.6 mm Rounded tips, 22 mm Curved 'Dolphin Nose' blades, 120 mm Long	MMSU1213S	MMSU1213CS
Westcott Scissors Straight Pointed	Pointed tips, 17 mm Straight blades, 120 mm Long	MMSU1214S	MMSU1214CS
Castroviejo Scissors Straight Rounded	0.5 mm x 0.6 mm Rounded tips, 11.5 mm Straight blades, 34 mm Handle	MMSU1217S	MMSU1217CS
Castroviejo Scissors Curved Rounded	0.5 mm x 0.6 mm Rounded tips, 11.5 mm Curved blades, 34 mm Handle	MMSU1218S	MMSU1218CS
Vannas Scissors Curved	Pointed tips on 8 mm curved blades, 1.5 mm when closed, 34 mm Handle	MMSU1140S	MMSU1140CS
Vannas Scissors Straight	Pointed tips on 8 mm straight blades, 1.5 mm when closed, 34 mm Handle	MMSU1141S	MMSU1141CS
Ong Capsule Scissors	Pointed tips, 9.5 mm Angled blades, 1.9 mm x 2.1 mm Thick at hinge, 34 mm Handle	MMSU1116S	MMSU1116CS
IOL Cutter (Osher)	0.2 mm x 0.5 mm Blunt tips, 9.5 mm Serrated blades, 66 mm Handle	MMSU1459S	MMSU1459CS

Dissection Scissors



**Stevens Tenotomy
Curved Rounded**




**Iris Scissors
Curved Pointed**



**Stevens Tenotomy
Straight Rounded**

Dissection Scissors	Description	SKU #
Stevens Tenotomy Curved Rounded	1.5 mm × 1.5 mm Rounded tips, Curved blades, 110 mm Long	MMSU0208S
Stevens Tenotomy Curved Pointed	0.5 mm × 0.5 mm Pointed tips, Curved blades, 110 mm Long	MMSU0209S
Stevens Tenotomy Straight Pointed	0.5 mm × 0.5 mm Pointed tips, Straight blades, 110 mm Long	MMSU0213S
Stevens Tenotomy Straight Rounded	1.5 mm × 1.5 mm Rounded tips, Straight blades, 110 mm Long	MMSU0214S
Enucleation Scissors	3 mm × 2.5 mm Rounded tips, 66 mm Blades with a 13 mm steep curve, 130 mm Long	MMSU0239S
Metzenbaum Scissors Curved	1.9 mm × 0.5 mm Rounded tips, 32.5 mm Curved blades, 120 mm Long	MMSU0243S
Iris Scissors Curved Pointed	0.1 mm × 0.7 mm Pointed tips, 25 mm Curved blades, 115 mm Long	MMSU0244S
Iris Scissors Straight Pointed	0.1 mm × 0.5 mm Pointed tips, 32.5 mm Straight blades, 120 mm Long	MMSU0249S
Iris Scissors Curved Rounded	1 mm × 1.5 mm Rounded tips, 30 mm Curved blades, 115 mm Long	MMSU0294S
Iris Scissors Straight Rounded	1 mm × 1.5 mm Rounded tips, 30 mm Straight blades, 115 mm Long	MMSU0299S

Dressing Scissors



Dressing Scissors



Dressing Scissors	Description	SKU #
Dressing Scissors Straight Pointed	0.3 mm Pointed tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle	MMSU0454S
Dressing Scissors Straight Rounded	0.8 mm Rounded tips, 24.5 mm Straight blades, 102 mm Long, Plastic handle	MMSU0455S
Dressing Scissors Curved Pointed	0.3 mm Pointed tips, 24.5 mm Curved blades, 100 mm Long, Plastic handle	MMSU0441S
Dressing Scissors Curved Rounded	0.7 mm Rounded tips, 24.5 mm Curved blades, 100 mm Long, Plastic handle	MMSU0440S

Speculums

A complete range of single-use speculums suitable for all ophthalmic procedures. Sprung-wire and adjustable-body designs with solid or open wire blades.



Adjustable Speculums

Adjustable Speculums	Description	SKU #
Mechanical Speculum Solid Blade	6 mm × 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1488S
Mechanical Speculum Open Blade	5 mm × 14 mm Open blades, 1 mm Wire, 38 mm Wide	MMSU1290S
Mechanical Speculum 45 Degree Solid Blade	7.5 mm × 14.5 mm 45 Degree angled polycarbonate solid blades, 0.8 mm Thick, Polycarbonate body and screw mechanical speculum	MMSU1566S
Mechanical Speculum V-Shaped Open Blade	5.5 mm × 18.5 mm Open v-shaped blades, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1494S
Adjustable Speculum	6.5 mm × 14 mm Open blades, 0.9 mm Wire, 25 mm Wide adjustable to 45 mm wide	MMSU1264S
Clarke Speculum	8.5 mm × 14.5 mm Open blades, 1.5 mm Wire, Metal body and screw mechanical speculum	MMSU1258S
Barrett Speculum	6 mm × 14.5 mm Curved blades, 32 mm Wide with screw locking mechanism	MMSU1266S



Barrett Speculum




Adjustable Speculum



Mechanical Speculum
Solid Blade

Sprung Wire Speculums




Kratz-Barraquer Speculum



Barraquer Speculum Flat Blade



Kratz-Barraquer Speculum Curved



Barraquer Speculum Solid Curved Blade



Barraquer Speculum Solid 18 mm Blade


Sprung Wire Speculums	Description	SKU #
Barraquer Speculum Solid Curved Blade	6 mm x 14.5 mm Solid curved polycarbonate blades, 0.7 mm Thick, 0.9 mm Wire, 30 mm	MMSU1474S
Barraquer Speculum Solid Straight Blade	7 mm x 14.5 mm Solid straight polycarbonate blades, 0.8 mm Thick, 0.9 mm Wire, 30 mm Wide	MMSU1475S
Barraquer Speculum Solid 18 mm Blade	6 mm x 14.5 mm Solid curved polycarbonate blades, 0.8 mm Thick, 1 mm Wire	MMSU1479S
Kratz-Barraquer Speculum	6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide	MMSU1252S
Kratz-Barraquer Speculum Temporal	6.5 mm x 14 mm Open blades, 0.9 mm Wire, 30 mm Wide	MMSU1253S
Kratz-Barraquer Speculum Curved	6.5 mm x 14.5 mm Open curved blades, 0.9 mm Wire, 30 mm Wide	MMSU1287S
Barraquer Speculum Flat Blade	6.5 mm x 14.5 mm Straight closed blades, 1 mm Wire, 30 mm Wide	MMSU1275S
Barraquer Speculum Temporal	6.5 mm x 14.5 mm Closed blades, 1 mm Wire, 30 mm Wide	MMSU1246S
Barraquer Speculum Curved Blade	5 mm x 16.5 mm Curved closed blades, 1 mm Wire, 38 mm Wide	MMSU1250S
Corcelle Speculum Solid Blade	6.5 mm x 8.5 mm Solid blades, 0.7 mm Thick, 55 mm Wide	MMSU1288S

Paediatric and Neonatal Speculums

Speculums Paediatric and Neonatal	Description	SKU #
Paediatric Mechanical Speculum Open	5 mm x 10.5 mm Open blades, 1 mm Wire, Polycarbonate body and screw mechanical speculum	MMSU1265S
Alphonso Paediatric Speculum	7 mm x 5.5 mm, 8 mm x 8.5 mm Closed blades, 1 mm Wire, 48 mm Wide	MMSU1260S
Neonatal Speculum Large	4 mm x 10 mm Closed blades, 0.8 mm Wire, 25 mm Wide	MMSU1261S
Neonatal Speculum Small	3 mm x 7 mm Closed blades, 0.8 mm Wire, 25 mm Wide	MMSU1262S
Paediatric Speculum	5.5 mm x 13.5 mm Closed blades, 0.9 mm Wire, 31.5 mm Wide	MMSU1278S
Cook Neonatal Speculum	5.5 mm x 8 mm Solid curved blades, 1 mm Thick, with Locking screw, Adjustable to 38 mm wide	MMSU1273S
Sauer Neonatal Speculum	4.6 mm x 6.6 mm Solid blades, 0.6 mm Wire spring arms, 18.5 mm Wide	MMSU1259S



Neonatal Speculum Small



Alphonso Paediatric Speculum



Cook Neonatal Speculum



Sauer Neonatal Speculum

Hooks and Retractors



Hooks and Retractors	Description	SKU #
Chevasse Muscle Hook 'W' Shaped	12 mm 'W' shaped hook at 90 degrees, 76 mm Handle	MMSU0232S
Curved Muscle Hook	10.3 mm Curved hook, 78 mm Handle	MMSU0253S
Desmarres Lid Retractor Small	11 mm × 11 mm × 5 mm Retractor tip, 35 mm Straight shaft, 76 mm Handle	MMSU3001S
Desmarres Lid Retractor Medium	14.5 mm × 10 mm × 8 mm Retractor tip, 30 mm Straight shaft, 76 mm Handle	MMSU3002S
Fisons Retractor	10.8 mm × 20 mm Flat retractor blade, 3 mm Lip, 76 mm Handle	MMSU3038S
Gillies Skin Hook	1 mm × 3 mm Hook, 31 mm Shaft, 76 mm Handle	MMSU3037S
Eye Guard	18 mm × 6 mm Rounded end, 20.5 mm × 6.5 mm Rounded end, 100 mm Long	MMSU1425S

Depressors



Depressors	Description	SKU #
Schocket Double-Ended Scleral Depressor	2.7 mm × 6.5 mm T shaped tip, 5.5 mm × 2.2 mm Rounded tip, 68 mm Handle	MMSU1405S
VR Double-Ended Depressor/Indentor	1.6 mm Indentor on 2 mm depressor, 1.6 mm Indentor on 3 mm depressor, 65 mm Handle	MMSU1407S
Curved Scleral Depressor	1.5 mm × 6.1 mm T shaped tip, 21 mm 50 Degree angled curved shaft, 78 mm Handle	MMSU1409S
Triangular Scleral Depressor	7 mm × 7 mm Triangular shaped tip, 0.75 mm Wire, 15 mm Straight shaft, 76 mm Handle	MMSU1247S
Ball Tip Scleral Depressor	3 mm Ball tip, 21.5 mm Straight shaft, 50 mm Handle	MMSU1408S

Lacrimal



Punctal Dilator



Lacrimal Probe



Golf Club Spud

Lacrimal	Description	SKU #
Punctum Seeker	5 mm Pointed tip, 18 mm Straight shaft, 77 mm Handle	MMSU1183S
Punctum Dilator/Seeker	0.45 mm x 21.5 mm Tip, 5 mm Pointed tip, 18 mm Straight shaft, 77 mm Handle	MMSU1193S
Punctum Dilator Small	0.3 mm x 20 mm Dilator, 30 mm Handle	MMSU1426S
Punctum Dilator Medium	0.55 mm x 24 mm Tip, 76 mm Handle	MMSU1180S
Punctum Dilator Fine/Medium	0.4 mm x 21.5 mm Tip, 0.6 mm x 18 mm Tip, 77 mm Handle	MMSU1181S
Lacrimal Probe	1 mm Pointed tip, 0.6 mm 80 mm Probe, 15 mm Handle	MMSU1427S
Punctum Dilator/Lacrimal Probe	0.28 mm x 25 mm Dilator, 0.7 mm x 44.5 mm Probe, 70 mm Handle	MMSU1419S
Bowmans Probes Size 000-00	0.7 mm (000G) x 50 mm Probe, 0.8 mm (00G) x 50 mm Probe, 30 mm Handle	MMSU1420S
Bowmans Probes Size 00-0	0.8 mm (00G) x 50 mm Probe, 0.9 mm (0G) x 50 mm Probe, 30 mm Handle	MMSU1421S
Golf Club Spud	0.1 mm x 0.8 mm x 2 mm Oval profile, 18 mm Straight shaft, 75 mm Handle	MMSU1417S




Markers and Gauges



Braunstein IVT Marker



Pointed IVT Marker



Double Ended IVT Guide Pad



Intravitreal Injection Guide

Markers and Gauges	Description	SKU #
Intravitreal Injection Guide	Polycarbonate lash guard with stainless-steel swage guide tube, Triangular base plate with 3 stabilizing studs, 39 mm, 45 Degree angled wishbone handle	MMSU1577S
Braunstein IVT Marker 3.5 4 mm Green	2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	MMSU1276S
Pointed IVT Marker 3.5 4 mm Blue	2 mm x 0.55 mm Tips with 3.5 mm spread, 2 mm x 0.65 mm Tips with 4 mm spread, Polycarbonate, 108 mm Long	MMSU1277S
Double Ended IVT Guide Pad 3.5 4 mm	7 mm Pad with 3.5 guide, 0.8 mm Pad with 4 mm guide, 6 mm 45 Degree angled shafts, 90 mm Long	MMSU1285S
Braunstein IVT Marker 3.5 4 mm	2 mm x 1.5 mm Tips with 3.5 mm spread, 2 mm x 1.5 mm Tips with 4 mm spread, 90 mm Long	MMSU1229S
Castroviejo Caliper Straight Tip	0.4 mm x 0.8 mm Straight tips, 30 mm Long, 0 mm to 25 mm in 0.5 mm Gradations	MMSU1401S
Castroviejo Caliper Angled Tip	0.25 mm x 0.3 mm Angled tips, 34 mm Long, 0 mm to 25 mm in 0.5 mm Gradations	MMSU1403S
Mendez Gauge	2.5 mm Polycarbonate ring with 13 mm ID, 0 through 180 Degree marks in increments of 10, 75 mm Handle	MMSU1152S
Toric Marker (Zuberbühler)	2.5 mm Wide polycarbonate ring with 13 mm ID, Marking blades at 0 and 180 degrees, Degree marks in increments of 10 degrees, 75 mm Handle	MMSU1245S
Horizontal Axis Marker	Spirit level guided horizontal marker, 75 mm Handle	MMSU1230S



Intravitreal Injection Guide

Clamps and Curettes



Clamps

Clamps	Description	SKU #
Chalazion Clamp Round	12 mm Round profile, 9 mm Round aperture, 90 mm Long	MMSU1150S
Chalazion Clamp Round Grooved	12 mm Round profile, 9 mm Round aperture, Grooved, 90 mm Long	MMSU1450S
Chalazion Clamp Oval	13 mm x 1.5 mm Rounded tips, Straight blades, 110 mm Long	MMSU1149S
Chalazion Clamp Oval Grooved	13 mm x 20 mm Oval profile, Grooved, 10 mm x 17 mm Aperture, 90 mm Long	MMSU1449S

Curettes



Curettes	Description	SKU #
Chalazion Curette 2 mm	2 mm Curette, 1.35 mm Depth, 115 mm Long	MMSU1155S
Chalazion Curette Double 1 mm 2 mm	1 mm Curette, 1 mm Depth, 2 mm Curette, 1.35 mm Depth, 128 mm Long	MMSU1157S
Chalazion Curette Double 2 mm 3 mm	2 mm Curette, 1.35 mm Depth, 3 mm Curette, 1.6 mm Depth, 128 mm Long	MMSU1158S

Flap-Lifters and Corneal Markers



LASIK Flap Elevator



LASIK Flap Separator



Y Shape
Corneal Marker



Flap-Lifters and Corneal Markers	Description	SKU #
LASIK Flap Elevator Curved 8 mm	0.2 mm × 0.15 mm Tip, Curved 8 mm tapered shaft at 90 degrees, 75 mm Handle	MMSU1177S
LASIK Flap Elevator Curved 10 mm	0.2 mm × 0.15 mm Tip, Curved 10 mm tapered shaft at 90 degrees, 75 mm Handle	MMSU1435S
LASIK Flap Separator and Elevator, Binder	0.1 mm × 0.3 mm Rounded tip 0.7 mm long, 10 mm Curved 60 degree angled shaft, 75 mm Handle	MMSU1206S
LASIK Flap Separator Flattened Tip	0.75 mm × 1.1 mm Flattened Tip, 23 mm Straight Shaft, 75 mm Handle	MMSU1176S
Double-Ring Corneal Marker	1.75 mm Ring overlapping 2 mm Ring, 2 mm Deep with ridged underside, 75 mm Handle	MMSU1189S
Y Shape Corneal Marker	2 mm × 8 mm Marking blades, 75 mm Handle	MMSU1134S
LASIK Flap Separator/Elevator (Seibel)	0.16 mm × 0.28 mm × 1 mm Tip, 20 mm Straight shaft, 0.5 mm × 10.5 mm pointed hook, 52 mm Handle	MMSU1172S
LASIK Flap Separator/Elevator	2.95 mm × 8 mm Triangular tip at a 60 degree angle, 24 mm Straight shaft, 0.53 mm Pointed tip, 9 mm 50 Degree angled shaft, 52 mm Handle	MMSU1205S
LASIK Flap Separator/Elevator Curved	0.15 mm × 0.3 mm × 8 mm Curved separator at 90 degrees to shaft, 0.15 mm × 0.3 mm × 1 mm Pocketing hook, 52 mm Handle	MMSU1294S
LASIK Flap Separator/Elevator Straight	0.15 mm × 0.3 mm × 8 mm Straight separator at 90 degree to shaft, 0.15 mm × 0.3 mm × 1 mm Pocketing hook, 52 mm Handle	MMSU1295S
LASIK Flap-Lifting Forceps	0.3 mm × 11.5 mm Tips, 45 Degree angled shaft, 80 mm Handle	MMSU1169S



LASIK Flap-Lifting Forceps

SMILE




SMILE Cornea	Description	SKU #
Reinstein Lenticule Separator	0.22 mm × 2 mm Pocketing hook, 1.2 mm × 2 mm Spatula tip, Curved separator at 90 degrees to shaft, 52 mm Handle	MMSU1297S
Reinstein Lenticule Separator (Hyperopic)	0.22 mm × 2 mm Pocketing hook, 1.2 mm × 2 mm Spatula tip, Curved separator at 90 degrees to shaft, 52 mm Handle	MMSU1497S
Pocketing Hook Separator	0.15 mm × 2 mm Pocketing hook angled at 90 degrees, 0.27 mm Tip increasing to 0.8 mm curved shaft separator, 8.5 mm Chord length, 90 Degrees to shaft, 52 mm Handle	MMSU1447S
IOL Removal Forceps	0.1 mm × 0.1 mm Ridges along 4 mm of 15 mm shaft, Angled at 30 degree, 52 mm Handle	MMSU1460S

Spatulas and Wells



Spatulas and Wells	Description	SKU #
Corneal Epithelium Spatula Flat	5 mm Spatula with 0.9 mm flat blade, 0.8 mm × 1.7 mm Round bladed tip, 75 mm Handle	MMSU1113S
Corneal Epithelium Spatula Round	2.65 mm × 5.5 mm × 0.35 mm Rounded blade, 75 mm Handle	MMSU1187S
Epithelium Hoe Flat	2.2 mm × 0.45 mm × 1.15 mm Micro hoe, 10.5 mm 45 Degree angled shaft, 75 mm Handle	MMSU1117S
Epithelium Hoe Curved	1.7 mm × 0.3 mm × 1 mm Micro hoe, 10.5 mm 45 Degree angled shaft, 75 mm Handle	MMSU1185S
Epithelium Well 6.5 mm	3 mm × 6.5 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1201S
Epithelium Well 7.5 mm	3 mm × 7.5 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1199S
Epithelium Well 8 mm	3 mm × 8 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1198S
Epithelium Well 8.5 mm	3 mm × 8.5 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1197S
Epithelium Well 9 mm	3 mm × 9 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1203S
Epithelium Well 9.5 mm	2 mm × 9.5 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1132S
Epithelium Well 9 mm Trepining	3 mm × 9 mm ID well, Angled 45 degrees, 75 mm Handle	MMSU1196S

Keratoplasty



Keratoplasty Instruments	Description	SKU #	SKU # CLIPS
Castroviejo Corneoscleral Scissors (L)	0.35 mm × 0.3 mm Pointed tips, 7 mm Curved blades, 34 mm Handle	MMSU1469S	MMSU1469CS
Castroviejo Corneoscleral Scissors (R)	0.35 mm × 0.3 mm Pointed tips, 7 mm Curved blades, 34 mm Handle	MMSU1470S	MMSU1470CS
Daya Descemets Scraper	0.6 mm × 0.6 mm Curved and reverse angled tip, 15 mm 45 Degree angled shaft, 75 mm Handle	MMSU1561S	-
'S' Marker	0.9 mm × 1.7 mm 'S' shaped marker, 0.2 mm Wide, 5 mm 45 Degree angled shaft, 75 mm Handle	MMSU1582S	-
Corneal Trephine 17 mm	17 mm Internal, 18.5 mm External, 20 mm High	MMSU0136S	-
Corneal Trephine 18 mm	18 mm Internal, 19 mm External, 30 mm High	MMSU0135S	-
0.1 mm Sinskey, 0.1 mm Reverse Sinskey	0.1 mm × 0.3 mm Tips, 11 mm 50 Degree angled shafts, 52 mm Handle	MMSU1466S	-

Corneal Ring




Corneal Ring	Description	SKU #
Corneal Ring Forceps	0.4 mm × 0.2 mm Triangular jaws, 10 mm Shaft angled at 35 degrees, 70 mm Handle	MMSU1462S
Corneal Ring Pocket Spreader	0.28 mm × 0.85 mm × 2.5 mm Oval profile tip, 18 mm Straight shaft, 75 mm Handle	MMSU1463S
Corneal Ring Pocketing Hook	0.2 mm × 0.9 mm Pointed tip, 18 mm Straight shaft, 75 mm Handle	MMSU1464S
Corneal Ring Pusher	0.55 mm × 3.5 mm Tip angled at 45 degrees, 18 mm Straight shaft, 75 mm Handle	MMSU1465S

Sclerectomy and Trabeculectomy



Sclerectomy and Trabeculectomy	Description	SKU #
Trabeculectomy Punch Titanium	1 mm x 0.9 mm Punch tip on 13 mm shaft, 88 mm Handle	MMSU1280S
Sclerectomy Forceps	0.2 mm x 1 mm Tips, 23 mm Straight shaft, 110 mm Long	MMSU1572S


I-Ring®



I-Ring	Description	SKU #
I-Ring Pupil Expander	Sterile, Single-use, 1/Ea	587001
I-Ring Pupil Expander	Sterile, Single-use, 6/Bx	587000
Sinsky Hook	Plastic handle, Angled 45 degrees, 10 mm from Tip, 0.85 mm Tip to bend distance, Overall length 130 mm	581442

Miscellaneous

Misc	SKU #
Clip - Long Handled Forceps (Pair)	MMSU1490S
Clip - Short Handled Forceps (Pair)	MMSU1491S
Clip - Sprung Scissors (Pair)	MMSU1492S
Clip - Micro Sprung Scissors (Pair)	MMSU1493S



Index

SKU	DESCRIPTION	PG	SKU	DESCRIPTION	PG
MMSU0068S	Straight Rampley	8	MMSU1139S	45° Toothed Forceps	5
MMSU0080S	Fine Stem Toothed Forceps	7	MMSU1140S	Vannas Scissors Curved	14
MMSU0085S	Epilation Forceps	7	MMSU1141S	Vannas Scissors Straight	14
MMSU0135S	Corneal Trephine 18 mm	24	MMSU1142S	Button Manipulator	12
MMSU0136S	Corneal Trephine 17 mm	24	MMSU1143S	McPherson Manipulator	13
MMSU0199S	Micro Cross-Action Bulldog Clamp	8	MMSU1144S	Drysdale Manipulator	12
MMSU0208S	Stevens Tenotomy Curved Rounded	15	MMSU1145S	45° Fine Toothed Forceps Short	5
MMSU0209S	Stevens Tenotomy Curved Pointed	15	MMSU1146S	Lens Loading Forceps	6
MMSU0213S	Stevens Tenotomy Straight Pointed	15	MMSU1147S	Sinskey Hook Straight	12
MMSU0214S	Stevens Tenotomy Straight Rounded	15	MMSU1149S	Chalazion Clamp Oval	21
MMSU0232S	Chevasse Muscle Hook 'W' Shaped	18	MMSU1150S	Chalazion Clamp Round Grooved	21
MMSU0238S	Adson Forceps, Toothed	7	MMSU1152S	Mendez Gauge	20
MMSU0239S	Enucleation Scissors	15	MMSU1153S	ICL Manipulator	13
MMSU0241S	Adson Forceps, Non-Toothed	7	MMSU1155S	Chalazion Curette 2 mm	21
MMSU0243S	Metzenbaum Scissors Curved	15	MMSU1157S	Chalazion Curette Double 1 mm 2 mm	21
MMSU0244S	Iris Scissors Curved Pointed	15	MMSU1158S	Chalazion Curette Double 2 mm 3 mm	21
MMSU0246S	Halstead Mosquito Straight	8	MMSU1168S	Toothed Forceps 15 mm Shaft	5
MMSU0247S	Iris Pointed Forceps	7	MMSU1169S	LASIK Flap-Lifting Forceps	22
MMSU0249S	Iris Scissors Straight Pointed	15	MMSU1172S	LASIK Flap Separator/Elevator (Seibel)	22
MMSU0253S	Curved Muscle Hook	18	MMSU1173S	Sinskey Hook 0.1 mm	12
MMSU0276S	Bishops Forceps	7	MMSU1176S	LASIK Flap Separator Flattened Tip	22
MMSU0277S	Moorfields Forceps	7	MMSU1177S	LASIK Flap Elevator Curved 8 mm	22
MMSU0278S	Iris Fine Toothed Forceps	7	MMSU1179S	Mushroom Manipulator	12
MMSU0279S	St Martins Forceps	7	MMSU1180S	Punctum Dilator Medium	19
MMSU0283S	Jewellers Forceps Fine Angled	7	MMSU1181S	Punctum Dilator Fine/Medium	19
MMSU0284S	Jewellers Forceps Curved/Plain	7	MMSU1183S	Punctum Seeker	19
MMSU0286S	Jewellers Forceps Fine Straight	7	MMSU1185S	Epithelium Hoe Curved	23
MMSU0294S	Iris Scissors Curved Rounded	15	MMSU1187S	Corneal Epithelium Spatula Round	23
MMSU0299S	Iris Scissors Straight Rounded	15	MMSU1189S	Double-Ring Corneal Marker	22
MMSU0440S	Dressing Scissors Curved Rounded	15	MMSU1193S	Punctum Dilator/Seeker	19
MMSU0441S	Dressing Scissors Curved Pointed	15	MMSU1196S	Epithelium Well 9 mm Trephining	23
MMSU0454S	Dressing Scissors Straight Pointed	15	MMSU1197S	Epithelium Well 8.5 mm	23
MMSU0455S	Dressing Scissors Straight Rounded	15	MMSU1198S	Epithelium Well 8 mm	23
MMSU1102S	Govan Manipulator	12	MMSU1199S	Epithelium Well 7.5 mm	23
MMSU1103S	Sinskey Hook 0.2 mm	12	MMSU1201S	Epithelium Well 6.5 mm	23
MMSU1104S	Sadiq Chopper	10	MMSU1203S	Epithelium Well 9 mm	23
MMSU1105S	Claw Chopper	10	MMSU1204S	Horizontal Y Manipulator	12
MMSU1106S	Rosen Chopper	10	MMSU1205S	LASIK Flap Separator/Elevator	22
MMSU1108S	Phaco Chopper	10	MMSU1206S	LASIK Flap Separator and Elevator, Binder	22
MMSU1109S	Vertical Y Manipulator	12	MMSU1207S	Lester Hook	13
MMSU1110S	Corneal Needle-Holder	9	MMSU1209S	Colibri Notched Forceps Short	5
MMSU1113S	Corneal Epithelium Spatula Flat	23	MMSU1210S	Curved Tying Forceps	6
MMSU1114S	Fine Notched Forceps Short (Hoskins)	5	MMSU1211S	Westcott Scissors Curved Pointed	14
MMSU1115S	Notched Forceps (Hoskins)	5	MMSU1212S	Westcott Scissors Straight Rounded	14
MMSU1116S	Ong Capsule Scissors	14	MMSU1213S	Westcott Scissors Curved Dolphin Nose	14
MMSU1117S	Epithelium Hoe Flat	23	MMSU1214S	Westcott Scissors Straight Pointed	14
MMSU1118S	Koch Chopper 1.9 mm	10	MMSU1217S	Castroviejo Scissors Straight Rounded	14
MMSU1119S	Greens Chopper Right Handed	10	MMSU1218S	Castroviejo Scissors Curved Rounded	14
MMSU1120S	Greens Chopper Left Handed	10	MMSU1219S	Straight Tying Forceps Short	6
MMSU1123S	Fine Toothed Forceps Short (Bonn)	5	MMSU1221S	Phaco Chopper Mushroom Manipulator	11
MMSU1124S	Toothed Forceps (Bonn)	5	MMSU1225S	Kuglen Hook Straight	13
MMSU1125S	Straight Tying Forceps	6	MMSU1228S	Koch Chopper 1.4 mm	10
MMSU1126S	IOL Implantation Forceps	6	MMSU1229S	Braunstein IVT Marker 3.5 4 mm	20
MMSU1127S	Angled Tying Forceps (McPherson)	6	MMSU1230S	Horizontal Axis Marker	20
MMSU1129S	Colibri Toothed Forceps Short	5	MMSU1231S	Horizontal Chopper/Kuglen Hook	11
MMSU1131S	Lens Folding Forceps Triangular	6	MMSU1238S	Ursell Chopper Manipulator	11
MMSU1132S	Epithelium Well 9.5 mm	23	MMSU1239S	Diablo Chopper Flat Button Manipulator	11
MMSU1134S	Y Shape Corneal Marker	22	MMSU1243S	Troutman Needle-Holder	9
MMSU1137S	Thornton Fixation Ring	7	MMSU1245S	Toric Marker (Zuberbühler)	20
MMSU1138S	Lens Folding Forceps Rounded	6	MMSU1246S	Barraquer Speculum Temporal	17

SKU	DESCRIPTION	PG	SKU	DESCRIPTION	PG
MMSU1247S	Triangular Scleral Depressor	18	MMSU1462S	Corneal Ring Forceps	24
MMSU1250S	Barraquer Speculum Curved Blade	17	MMSU1463S	Corneal Ring Pocket Spreader	24
MMSU1252S	Kratz-Barraquer Speculum	17	MMSU1464S	Corneal Ring Pocketing Hook	24
MMSU1253S	Kratz-Barraquer Speculum Temporal	17	MMSU1465S	Corneal Ring Pusher	24
MMSU1258S	Clarke Speculum	16	MMSU1466S	0.1 mm Sinskey, 0.1 mm Reverse Sinskey	24
MMSU1259S	Sauer Neonatal Speculum	17	MMSU1469S	Castroviejo Corneoscleral Scissors (L)	24
MMSU1260S	Alphonso Paediatric Speculum	17	MMSU1470S	Castroviejo Corneoscleral Scissors (R)	24
MMSU1261S	Neonatal Speculum Large	17	MMSU1471S	Denis Nucleus Chopper	10
MMSU1262S	Neonatal Speculum Small	17	MMSU1472S	Jewellers Forceps Angled/Smooth	7
MMSU1264S	Adjustable Speculum	16	MMSU1473S	0.1 mm Reverse Sinskey	13
MMSU1265S	Paediatric Mechanical Speculum Open	17	MMSU1474S	Barraquer Speculum Solid Curved Blade	17
MMSU1266S	Barrett Speculum	16	MMSU1475S	Barraquer Speculum Solid Straight Blade	17
MMSU1273S	Cook Neonatal Speculum	17	MMSU1479S	Barraquer Speculum Solid 18 mm Blade	17
MMSU1275S	Barraquer Speculum Flat Blade	17	MMSU1480S	Reverse Dialling Hook	13
MMSU1276S	Braunstein IVT Marker 3.5 4 mm Green	20	MMSU1487S	Phaco Chopper Mushroom Manipulator Fine	11
MMSU1277S	Pointed IVT Marker 3.5 4 mm Blue	20	MMSU1488S	Mechanical Speculum Solid Blade	16
MMSU1278S	Paediatric Speculum	17	MMSU1489S	Neuhann Nucleus Divider	10
MMSU1280S	Trabeculectomy Punch Titanium	25	MMSU1490S	Clip - Long Handled Forceps (Pair)	25
MMSU1284S	Contact Lens Forceps	7	MMSU1491S	Clip - Short Handled Forceps (Pair)	25
MMSU1285S	Double Ended IVT Guide Pad 3.5 4 mm	20	MMSU1492S	Clip - Sprung Scissors (Pair)	25
MMSU1287S	Kratz-Barraquer Speculum Curved	17	MMSU1493S	Clip - Micro Sprung Scissors (Pair)	25
MMSU1288S	Corcelle Speculum Solid Blade	17	MMSU1494S	Mechanical Speculum V-Shaped Open Blade	16
MMSU1290S	Mechanical Speculum Open Blade	16	MMSU1495S	Fixation Ring	7
MMSU1291S	Spatula	12	MMSU1496S	Phaco Chopper Govan Manipulator	11
MMSU1294S	LASIK Flap Separator/Elevator Curved	22	MMSU1497S	Reinstein Lenticule Separator (Hyperopic)	23
MMSU1295S	LASIK Flap Separator/Elevator Straight	22	MMSU1561S	Daya Descemet's Scraper	24
MMSU1297S	Reinstein Lenticule Separator	23	MMSU1566S	Mechanical Speculum 45 Degree Solid Blade	16
MMSU1401S	Castroviejo Caliper Straight Tip	20	MMSU1567S	MICS Capsulorhexis - Curved/Sharp	4
MMSU1403S	Castroviejo Caliper Angled Tip	20	MMSU1568S	MICS Capsulorhexis - Straight/Sharp	4
MMSU1405S	Schocket Double-Ended Scleral Depressor	18	MMSU1571S	MICS Capsulorhexis - Straight/Grasping	4
MMSU1407S	VR Double-Ended Depressor/Indentor	18	MMSU1572S	Sclerectomy Forceps	25
MMSU1408S	Ball Tip Scleral Depressor	18	MMSU1573S	MICS Capsulorhexis - Curved/Grasping	4
MMSU1409S	Curved Scleral Depressor	18	MMSU1577S	Intravitreal Injection Guide	20
MMSU1411S	Scleral Plug Forceps	7	MMSU1578S	Ramley Short	8
MMSU1414S	Birks Tying Forceps	6	MMSU1582S	'S' Marker	24
MMSU1417S	Golf Club Spud	19	MMSU1583S	Waqar Suture Removal Forceps	6
MMSU1419S	Punctum Dilator/Lacrimal Probe	19	MMSU1592S	Barraquer Iris Spatula Round	13
MMSU1420S	Bowmans Probes Size 000-00	19	MMSU3001S	Desmarres Lid Retractor Small	18
MMSU1421S	Bowmans Probes Size 00-0	19	MMSU3002S	Desmarres Lid Retractor Medium	18
MMSU1425S	Eye Guard	18	MMSU3007S	Halstead Mosquito Fine Curved	8
MMSU1426S	Punctum Dilator Small	19	MMSU3008S	Halstead Mosquito Fine Straight	8
MMSU1427S	Lacrimal Probe	19	MMSU3009S	Castroviejo Needle-Holder Straight	9
MMSU1428S	Diablo Manipulator	12	MMSU3011S	Castroviejo Needle-Holder Curved	9
MMSU1431S	Koch Chopper 1.4 mm Button Manipulator	11	MMSU3012S	Castroviejo Needle-Holder Straight with Lock	9
MMSU1432S	Flat Round Tip Manipulator	12	MMSU3013S	Castroviejo Needle-Holder Curved with Lock	9
MMSU1433S	Kuglen Hook Angled	13	MMSU3015S	Barraquer Needle-Holder	9
MMSU1434S	Butcher Rotator	12	MMSU3017S	Corneal Needle-Holder Long	9
MMSU1435S	LASIK Flap Elevator Curved 10 mm	22	MMSU3037S	Gillies Skin Hook	18
MMSU1436S	Hagège Manipulator	12	MMSU3038S	Fisons Retractor	18
MMSU1440S	Koch Chopper 1.4 mm Thin Shaft	10	581441	Micro Manipulator	13
MMSU1441S	Nagahara Chopper	10	581442	Sinskey Hook	13, 25
MMSU1445S	Max Fine Micro-Corneal Forceps	7	581486	Next Generation Capsulorhexis Forceps Grasping	4
MMSU1447S	Pocketing Hook Separator	23	581488	Next Generation Capsulorhexis Forceps Sharp	4
MMSU1449S	Chalazion Clamp Oval	21	587000	I-Ring Pupil Expander 6/Bx	25
MMSU1450S	Chalazion Clamp Round Grooved	21	587001	I-Ring Pupil Expander 1/Ea	25
MMSU1452S	Ayliffe Chopper Chopper	11			
MMSU1457S	Barraquer Iris Spatula	12			
MMSU1458S	Lebuisson Manipulator	13			
MMSU1459S	IOL Cutter (Osher)	14			
MMSU1460S	IOL Removal Forceps	6, 23			



Beaver-Visitec Sales International LTD

Tel: + 44 1865 601 256 Option 3
Fax: + 44 1865 595 761
Email: ordersuk@bvimedical.com

Belgium

Tel: + 32 52 280518 Option 3
Fax: + 32 52 890243
Email: bestellingen@bvimedical.com

France, Luxembourg

Tel: + 33 4 85880033 Option 3
Fax: + 33 4 88714057
Email: commandes@bvimedical.com

Germany, Switzerland, Austria

Tel: + 49 622 16479854 Option 3
Fax: + 49 622 13929883
Email: bestellungen@bvimedical.com

Italy

Tel: + 39 02 89877520 Option 3
Fax: + 39 02 44386868
Email: ordini@bvimedical.com

Japan

Tel: + 03 5651 2860
Fax: + 03 3662 2242
Email: bvi.cs@bvimedical.com

Netherlands

Tel: + 31 20808 3202 Option 3
Fax: + 31 20203 1293
Email: bestellingen@bvimedical.com

Spain

Tel: + 34 91 1230335 Option 3
Fax: + 34 91 1849957
Email: ordenes@bvimedical.com

All Other Regions

Tel: + 44 1865 601.256 Option 3
Fax: + 44 1865 595 761
Email: ordersrow@bvimedical.com



bvimedical.com



BVI
bvimedical.com